

What is it for?

Plexr can treat a variety of skin imperfections, by performing treatments for the removal of benign neoformations and elastic retraction treatments, thus generating skin lifting. Plexr is the gold standard in non-invasive surgery.

How does it work?

Plexr is a device that ionizes the air present between its terminal tip and the epidermis, thus generating a micro-plasma beam. In contact with the epidermis the micro-plasma beam sublimates the corneocytes, i.e. the upper cells that form the protective layer of the epidermis. At the same time, an amount of energy selectively reaches the right depth of the dermis, without irradiation (laser/light) or electric shock (electrosurgical/current). This stimulation allows an immediate contraction and reorganization of the collagen fibres, a formation of new collagen and a renewal of the epidermal tissue.

How long is the treatment?

Depending on the type of procedure – treatment time ranges between 10-60 minutes.

Are there any contraindications?

There are no permanent side effects for treatments performed with Plexr. Scabs, swelling and redness are part of the normal healing process. This normally lasts 7-20 days.

How long does the results last?

Results are immediately visible and absolutely lasting over time. For some types of treatment it may take more than one session to maximise the effectiveness of the results.

TREATMENT

PRICE

Upper & Lower Blepharoplasty	R4400
Upper OR Lower Blepharoplasty	R2200
Peri Oral (Smokers Lines)	R2400
Ear Wrinkles	R2400
Mini Facelift	R7000
Active Acne and Acne Scarring	From R950
Scarring	From R750
Skin Spots (Warts, Freckles, Moles, Xanthelasma)	From R750
Tattoo Removal	From R1200
Gummy Smile (at dental practice)	R1500
Tongue/Lip Tie (at dental practice)	R1500
Shower Facial (includes meso.pro solution)	R1200